

Germline Testing for Hereditary Breast and Gynecologic Cancer

REFERENCE GUIDE

Genetic Testing For Hereditary Breast and Gynecologic Cancer is Recommended By:

American Society of Breast Surgeons (ASBrS)¹

American College of Obstetricians and Gynecologists (ACOG)²

Society of Gynecologic Oncology (SGO)²

NCCN Clinical Practice Guidelines In Oncology (NCCN Guidelines[®])³

GUIDELINES RECOMMEND GENETIC TESTING IF YOUR PATIENT HAS A PERSONAL OR FAMILY HISTORY OF ANY OF THE FOLLOWING SIGNS* OF HEREDITARY CANCER:

- ✓ **MULTIPLE CANCERS**
 - › 2 or more primary cancers in the same person
 - › 3 or more cancers on the same side of the family

- ✓ **EARLY-ONSET CANCERS**
 - › Breast cancer diagnosed before 45
 - › Triple negative breast cancer before 60
 - › Uterine cancer diagnosed before 50

- ✓ **RARE CANCERS**
 - › Male breast
 - › Ovarian
 - › Pancreatic
 - › Metastatic prostate

- ✓ **ABNORMAL TUMOR SCREENING**
 - › Tumors with microsatellite instability (MSI) or loss of immunohistochemical (IHC) staining

- ✓ **ANCESTRY**
 - › Ashkenazi Jewish with breast cancer

*Adapted from national published guidelines

PALB2 MUTATION

Case Example

Screening and surgical implications
for patient & family members



IMPACT ON PATIENT MEDICAL MANAGEMENT

- › Begin annual mammograms and breast MRIs³ (starting now instead of at 40)
- › Consider risk-reducing mastectomy depending on family history³
- › Consider risk-reducing salpingo-oophorectomy based on family history³



IMPACT ON FAMILY MEMBERS

- › Genetic testing recommended for family members³
- › If positive, medical management of family members may be impacted. For example:
 - › Possible increased risk for second primary breast cancer in mother (may consider bilateral mastectomy)
 - › Mother may consider risk-reducing salpingo-oophorectomy



31 year-old female

Personal History

No personal history of cancer

Family History

Mother

Breast Cancer (55)

Maternal Aunt

Ovarian Cancer (70)

Paternal Grandmother

Breast Cancer (68)



Patient meets *BRCA1/2* genetic testing criteria




Ordered **BRCANext™**



Genetic testing identified a **PALB2** mutation




Results Interpretation and Counseling to Inform Medical Management



GENETIC TESTING RESULTS MAY ALSO INFORM TREATMENT RECOMMENDATIONS

- › PARP inhibitor therapy for *BRCA1* or *BRCA2* positive patients with ovarian cancer⁴⁻¹²
- › Immunotherapy for patients with a germline mutation in a mismatch repair gene¹³

HEREDITARY BREAST AND GYNECOLOGIC CANCER TESTING OPTIONS



CancerNext®

36-gene test covering a range of cancers; NCCN Guidelines® provide recommendations regarding risk management for most genes in the panel³




BRCANext-Expanded™

23-gene test for hereditary breast and gynecologic cancers; NCCN Guidelines provide recommendations regarding risk management for most genes in the panel³



BRCANext™

18-gene test for hereditary breast and gynecologic cancers; NCCN Guidelines provide recommendations regarding risk management for **all genes** in the panel³



BRCAplus®

8-gene STAT, breast cancer test; NCCN Guidelines provide recommendations regarding risk management for **all genes** in the panel³



Add +RNAinsight® to Find More Mutations, Decrease Variants of Unknown Significance, and Provide More Accurate Results¹⁴⁻¹⁶

+RNAinsight is Available With Any Hereditary Cancer Panel*

➤ SAMPLE REQUIREMENTS

One kit that includes 1 EDTA tube (DNA) and 1 PAXgene® tube (RNA)

➤ NO IMPACT ON TURNAROUND TIME OR COST

DID YOU KNOW?

Adding RNA genetic testing decreases variants of unknown significance by **5%** and increases diagnostic yield by **7%**¹⁶

*Exceptions: BRCAplus® and TumorNext® tests

References

1. Manahan ER, et al. *Ann Surg Oncol*. 2019 Oct;26(10):3025-3031. 2. *Obstetrics & Gynecology*. 2017 Sep;130(3):e110-e126. 3. Referenced with permission from the NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines®) for Genetic/Familial High-Risk Assessment: Breast, Ovarian, and Pancreatic. V1.2020. © National Comprehensive Cancer Network, Inc. 2019. All rights reserved. Accessed June 11, 2020. To view the most recent and complete version of the guideline, go online to NCCN.org. NCCN makes no warranties of any kind whatsoever regarding their content, use or application and disclaims any responsibility for their application or use in any way. 4. Banerjee S & Kaye S. *Curr Oncol Rep*. 2011 Dec;13(6):442-9. 5. Burgess M & Puhalla S. *Front Oncol*. 2014 Feb 27;4:19. 6. Yamamoto KN, et al. *PLoS One*. 2014 Aug 26;9(8): e105724. 7. Moore, et al. *NEJM* 2018 Dec 27; 379(26):2495-505. 8. Ledermann, et al. *Lancet Oncol*. 2014;15(8):852-861. 9. Pujade-Lauraine, et al. *Lancet Oncol*. 2017;18:1274-1284. 10. Mirza, et al. *N Engl J Med*. 2016;375:2154-2164. 11. Coleman RL, et al. *Lancet*. 2017 Oct 28;390(10106):1949-1961. 12. Swisher, et al. *Lancet Oncology* 2017;18: 75-87 13. Zhao P, Li, et al. *J Hematol Oncol* 2019 May;12:54. 14. Landrith T et al. *npj Precision Oncology*. 2020. 15. Karam et al. *JAMA Network Open*. 2019. 16. Ambry Genetics, internal data on file from first 2,500 cases tested with +RNAinsight.

BREAST AND GYNECOLOGIC GENES AND ASSOCIATED RISKS

GENE(S)	ASSOCIATED CANCERS AND RISKS
<i>ATM</i> *	Breast (2-4 fold), pancreatic, prostate
<i>BRCA1</i> *	Female breast (57-87%), ovarian (39-40%), pancreatic, melanoma, prostate, male breast
<i>BRCA2</i> *	Female breast (45-84%), ovarian (11-18%), pancreatic, melanoma, prostate (15%), male breast (>6%)
<i>CDH1</i> *	Diffuse gastric (67-83%), female lobular breast (39-52%)
<i>CHEK2</i> *	Breast (2 fold), colorectal, prostate, other
<i>PALB2</i> *	Breast (33-58%), pancreatic, ovarian, possibly prostate, male breast
<i>PTEN</i> *	Breast (25-85%), thyroid (10-35%), uterine (5-28%), melanoma (up to 6%), other
<i>TP53</i> *	Breast, sarcoma, brain, adrenocortical, leukemia, other
<i>BRIP1</i> *	Ovarian (up to 9%), breast
<i>EPCAM</i>	Colorectal (52-82%), uterine (12-55%), possibly prostate, other
<i>MLH1</i>	Colorectal (52-82%), uterine (25-60%), stomach (6-13%), ovarian (4-12%), prostate (2 fold), other
<i>MSH2</i>	Colorectal (52-82%), uterine (25-60%), stomach (6-13%), ovarian (4-12%), prostate (2 fold), other
<i>MSH6</i>	Colorectal (20-44%), uterine (up to 44%), prostate (2 fold), other
<i>NBN</i> *	Breast, possibly ovarian, prostate
<i>NF1</i> *	Female breast (3-5 fold), malignant nerve sheath tumors (8-13%), PGL/PCC (up to 7%), brain, other
<i>PMS2</i> *	Colorectal (15-20%), uterine (15%), possibly prostate, other
<i>RAD51C</i> *	Ovarian (5-9%), breast
<i>RAD51D</i> *	Ovarian (10-12%), breast, prostate
<i>BARD1</i>	Breast, possibly ovarian
<i>DICER1</i>	Pleuropulmonary blastoma, cystic nephroma, ovarian sex-cord tumors, brain, other
<i>RECQL</i>	Breast
<i>SMARCA4</i>	Ovarian (small cell carcinoma, hypercalcemic type), brain, other
<i>STK11</i> *	GI cancers (up to 57%) breast (up to 45%), pancreatic, other
<i>APC</i> *	Colorectal (up to 100%), small bowel, stomach, pancreatic, thyroid, other
<i>AXIN2</i> *	Colorectal
<i>BMPRI1</i> */ <i>SMAD4</i> *	Colorectal (40-50%), stomach (up to 21%)
<i>CDK4</i>	Melanoma (up to 74% by age 50)
<i>CDKN2A</i> *	Melanoma (28-67%), pancreatic (17-25%)
<i>GREM1</i>	Colorectal
<i>HOXB13</i>	Prostate
<i>MSH3</i> *	Colorectal
<i>MUTYH</i> *	Biallelic mutations: colorectal (up to 80%), stomach, duodenal, uterine, Monoallelic mutations: colorectal (up to 2 fold)
<i>NTHL1</i> *	Colorectal
<i>POLD1</i> */ <i>POLE</i> *	Colorectal

BRCAPlus
BRCANext
BRCANext-Expanded
CancerNext

* Risk management recommendations are provided in the NCCN Guidelines³

+RNAinsight performed for this gene (RNA analysis available only for exons 1-10 for *PMS2*)